

APPENDIX J

10 Year Useful Life

- J.1.1 The Useful Life of the Andersen AN74i and AN74ix sterilizers is ten (10) years from the date of manufacture. The Lot/Serial Number on the back of the sterilizer indicates the year the sterilizer was manufactured. For example, the first two numbers of the Lot/Serial Number 283508 indicate the sterilizer was manufactured in the year 2008. Useful Life refers to the time period during which Andersen Sterilizers will maintain a spare parts inventory and provide services to repair and support the sterilizer. Useful Life also refers to the time period, after which the sterilizer should be removed from service and replaced with a newer sterilizer.
- J.1.2 If your sterilizer is in need of repair after the Limited One Year Warranty Period has expired and prior to the expiration of the Useful Life, you should contact your Andersen Customer representative at 800-523-1276 if you are located in the United States or Canada. If you are located outside of these two countries, you should contact your local distributor.